

Stryker
Power-PRO 2



Manual Stretcher & Power-LOAD


Table of Contents

- 1) Manual Stretcher (pg. 3)
- 2) Loading Manual Stretcher / Working Power-LOAD (pg. 5)
- 3) Unloading Manual Stretcher / Working Power-LOAD (pg. 10)
- 4) Loading / Power Loss of Power-LOAD (pg. 14)
- 5) Unloading / Power Loss of Power-LOAD (pg. 18)
- 6) Releasing the Power-LOAD Trolley (pg. 22)

Manual Stretcher

Manual Stretcher

1. Two operators are ready to support the weight at each end of the stretcher
2. Lift the weight off the stretcher wheels
3. Pull the red manual release handle at the foot-end and release to lock in position



Loading Manual Stretcher / Working Power-LOAD

Loading Manual Stretcher / Working Power-LOAD

1. Manually lift the stretcher to loading height
2. Lock the stretcher into the Power-LOAD – confirm the LED antler lights are green



Loading Manual Stretcher / Working Power-LOAD

3. Hold the UP arrow on the Power-LOAD control panel until the lifting arms stop



Loading Manual Stretcher / Working Power-LOAD

4. Pull the red manual release handle on the stretcher and manually retract the stretcher legs – **continue to squeeze the handle while you load the stretcher**



Loading Manual Stretcher / Working Power-LOAD

5. Pull the stretcher side-to-side to confirm it is locked in the track
6. If there is space between the track and stretcher (area circled below), pull the manual release handle once more



Unloading Manual Stretcher / Working Power-LOAD

Unloading Manual Stretcher / Working Power-LOAD

1. Unload the stretcher – confirm the Power-LOAD is locked and the LED antler lights are green
2. Pull the red manual release handle and manually extend the stretcher legs
3. Release the handle to lock the stretcher legs into position

Unloading Manual Stretcher / Working Power-LOAD

4. Once legs are fully extended, hold the DOWN arrow on the Power-LOAD control panel until the lifting arms stop



Unloading Manual Stretcher / Working Power-LOAD

5. Pull the red release handle at the head-end of the trolley to release the stretcher from the Power-LOAD



Loading Stretcher / Power Loss of Power-LOAD

Loading Stretcher / Power Loss of Power-LOAD

1. Lift the lifting arms and push the Power-LOAD trolley to the end of the track until locked into place



Loading Stretcher / Power Loss of Power-LOAD

2. Raise the stretcher to loading height
3. Hook the safety bar onto the safety hook at the head-end of the stretcher



Loading Stretcher / Power Loss of Power-LOAD

4. Lift the stretcher, retract the legs, and load into the ambulance
5. Guide the stretcher along the track and lock into the antlers of the trolley
6. Pull the foot-end of the stretcher side-to-side to ensure it is locked into place (LED antler lights will not work on the Power-LOAD)

Unloading Stretcher / Power Loss of Power-LOAD

Unloading Stretcher / Power Loss of Power-LOAD

Important!

Be ready to support the weight of the stretcher

The Power-LOAD will **not** be attached to the stretcher and will not assist with unloading

Unloading Stretcher / Power Loss of Power-LOAD

1. Unload the stretcher – go slow and ensure the safety bar hooks onto the yellow safety hook



Unloading Stretcher / Power Loss of Power-LOAD

2. Extend the stretcher legs to the floor
3. Release the safety bar from the hook by pulling the head-end release handle on the stretcher



Releasing the Power-LOAD Trolley

(If locked at the end of the track)

Releasing the Power-LOAD Trolley

1. Press and hold the release button located at the end of the trolley across from the jump seat
2. While holding the release button, lift the arms of the Power-LOAD and pull it out to release



Releasing the Power-LOAD Trolley

If the trolley does not release when trying to unload the stretcher, follow the same steps:

1. One partner presses and holds the trolley release button (across from the jump seat)
2. At the same time, the other partner begins to unload the stretcher.

Thank you!

If you have any questions, please reach out to your Education Team.

